

Instructions for use - Please read carefully!

# Lysoformin<sup>®</sup> rapid wipes

Pre-impregnated tissues for a rapid contact disinfection and cleaning of medical devices and disinfection of patient-near (small) surfaces in the medical and veterinary field

## Lysoformin<sup>®</sup> rapid wipes

Ready-to-use impregnated tissues in white stand-up plastic pouches

### Composition

100 tissues contain :

1,5 litre soaking solution of Lysoformin<sup>®</sup> rapid at 1%

Active substances per tissue :

8,85 mg g Didecyl dimethyl ammonium chloride

11,1 mg Bis-(3-aminopropyl) dodecylamine

### PROPERTIES

- 100 tissues (30 x 30 cm) per roll, one tissue is enough for a surface of about 1 m<sup>2</sup>
- economical in use, individual tear-off
- absorbent and lint-free fleece
- handy stand-up pouch with carrying handle, easy to use and transport
- EN 16615 (4-field test), product VAH listed

Surface disinfection and cleaning	
Prevention of nosocomial infections	5 mn contact time
Efficacy spectrum Virucidal against enveloped viruses	10 mn contact time

### RANGE OF APPLICATION

Rapid disinfection and cleaning of medical devices and patient-near (small) surfaces in the medical and veterinary areas.

*For professional use*

### MODE OF APPLICATION

Completely wet the surface to be disinfected and wait until the end of the prescribed contact time!

The user has to write the date of the first opening on the round label. Die impregnated tissues can be used up to 28 days after first opening the pouch.

#### Picture 1

Write date of the first opening on the round label, peel it off and stick it on the cap.

#### Picture 2

Put on protective gloves

#### Picture 3

Carefully cut out the bag foil inside the black screw ring with a sharp, pointed knife and take it out completely.

#### Picture 4

Pull out the first tissue though the removal cap completely and throw it away.

#### Picture 5

Close the lid of the removal cap and screw it onto the black thread.

### Preparation of the removal cap in case of re-use

An insufficient hygiene can lead to the contamination of the extraction system. Therefore the used extraction closure system should be cleaned every time the pouch is changed.

To do so unscrew the cap, open it and thoroughly clean it under running water. Pay particular attention to the thread. Then dry the removal cap with disposable paper towels and keep dry and protected until re-use. Alternatively, the removal cap can also be machine-processed (max 60 °C).

Before screwing the cap onto a new pouch it is recommended to disinfect the closure system with an alcoholic surface disinfectant.

### **WARNINGS / STORAGE**

- Use Lysoformin<sup>®</sup> rapid wipes carefully!  
Always read label and product information before use.
- Store Lysoformin<sup>®</sup> rapid wipes in the closed original container and in the upright position. Keep the product in a dry, cool but frost free and well ventilated place and out of reach of children.
- Wear adequate protective gloves when handling the produc

### **STANDARD PACKAGING / DELIVERY UNITS**

4 standard pouches containing 1 roll with 100 tissues each  
1 extracting closure system (cap)  
1 mini cutter

### **article numbers**

9235  
3405

**Lysoform Dr. Hans Rosemann GmbH**  
**Kaiser-Wilhelm-Str. 133, D – 12247 Berlin**  
**Tel.: +49-30-77992-0   Telefax: +49-30-77992-219**  
**[www.lysoform.de](http://www.lysoform.de), e-mail: [kontakt@lysoform.de](mailto:kontakt@lysoform.de)**

**2019/06/21\_02**