

Instructions for use - Please read carefully!

## Lysoformin® 3000

### Concentrate

For professional use

**Disinfecting and cleaning concentrate for instruments, equipment and surfaces; particularly appropriate for flexible endoscopes**

### COMPOSITION

100 g **solution** contains:

7.5 g Glyoxal, 9.5 g Glutaral, 9.6 g Didecyl Dimethyl Ammonium Chloride, Isotridecanol (ethoxylated)

### RANGE OF APPLICATION

- **Disinfection and cleaning of instruments, endoscopes**, as well as equipment and surfaces and textiles as well  
as medical devices) **in one Step only** in case of all bacterial (tb incl.) and viral infections
- Use according to the preparation instructions of the manufacturer of the medical device
- Control of nosocomial infections

### MODE OF APPLICATION / DOSAGE

**Lysoformin® 3000** is a liquid concentrate with which any use-dilution required can be made by simply adding water (at max. room temperature).

**Important:** Always fill container with tap water first and then add **Lysoformin® 3000**. Prepare use solution shortly before use (acc. to the Directives for hygiene at the hospital and prevention of infections, RKI). Solutions ready for use without organic soil (disinfection of pre-cleaned instruments) are efficacious at least 14 days

Disinfection of the medical instruments (incl. tb) by immersion / soaking		
Strong contamination including M. terrae	2.50 %	30 mn
	1.50 %	1 h
Virucidal against enveloped viruses* including HBV/HCV/HIV	0.25 %	15 mn
Virucidal*	2.00 %	30 mn
Polyoma virus (SV 40)	0.50 %	30 mn
Adeno virus	0.25 %	15 mn
Rota virus	0.25 %	1 mn
Noro virus (MNV)	1.00 %	5 mn
	0.50 %	15 mn
Textiles		
Mops, floor cloths	0.25 %	4 h
Scrubbing/wiping disinfection and cleaning of surfaces		
Prophylactic treatment of nosocomial infections	0.25 %	4 h
and routine disinfection	0.50%	1 h
	1.00%	15 mn

\*according to the recommendation of RKI 2017

**Important:** Place objects to be disinfected in the solution in such a way that all surfaces are in contact with the product. Close container tightly with lid. Rinse instruments with water after disinfection/cleaning. Use germ-free water to rinse endoscopes.

## Lysoform – Hygiene with System

Dosing system for 8 l of water \ Concentration	0.25 %	0.5%	1.0%	1.5 %	2.0 %
Dosing-bag (1 bag = 20 ml)	1 bag	2 bags	4 bags	6 bags	8 bags
Dosing pump (1 stroke = 20 ml)	1 stroke	2 strokes	4 strokes	6 strokes	8 strokes
<b>Lysoform dosing devices</b>	Adjust concentration once only to obtain the concentration wished for				

The necessary amount of water is calculated by subtracting the amount of disinfectant needed to obtain the concentration wished for from the total amount of solution needed

### **WARNINGS / STORAGE**

- Keep disinfectants in a cool but frost-free, well ventilated and dry place and out of reach of children
- Wear suitable protective clothing, adequate protective gloves and eye protection when handling the concentrated product.
- Use only in well-ventilated rooms
- Do not use beyond expiry date.
- Shelf life of the solution ready for use and free of organic residue\*: 14 days when the product is stored under appropriate conditions (in a cool and dry place) in clean and tightly closed containers.

### **PROPERTIES**

- Control of bacteria (tb incl.), fungus, viruses
- Good compatibility with **instruments including flexible endoscopes**, equipment and surfaces
- Fresh non-clinging scent
- Contains no formaldehyde, no phenols, no phosphates, no oxidizing compounds
- Cleans without adding detergents
- Resists hard water, biodegradable

### **SPECIAL PRECAUTIONS**

a) Incompatibilities:

Adding the usual detergents as well as rests of care products present on the instruments can impair or neutralize the disinfecting properties of **Lysoformin 3000®** or its compatibility with the material of the instruments.

b) Counter-measures in case of improper use:

Should **Lysoformin 3000®** come into contact with the mucous membranes (eyes, throat, stomach etc.) get medical attention immediately. As an emergency measure dilute product immediately: if it has been swallowed, drink **plenty** of water; after contact with the eyes, rinse thoroughly with **plenty** of water.

c) Use **Lysoformin 3000®** safely. Always read the label and product information before use.

Certified and mentioned in the relevant lists.

\* period of time during which the proven efficacy of the disinfectant solution stays intact, if the product was prepared in the concentration wished for according to the recommendations of the manufacturer, was then covered and left unused during the period of time recommended by the manufacturer.

### **STANDARD PACKAGING / TRANSPORT UNITS**

200 x 20 ml dosing bags = 4 litres in the cardboard box  
12 x 1000 ml bottles  
2 x 5 l-canisters  
1 x 20 l-canister

### **Article No:**

9508  
9256  
9265  
9330

**Accessories:**Dosing pump for 5 l can (1 stroke  $\cong$  20 ml)

3135

**Hygiene Technique:**

Lysoform Dosing Device

9174

**Lysoform Dr. Hans Rosemann GmbH**

Kaiser-Wilhelm-Str. 133, D-12247 Berlin

Tel.: +49-30-77992-0    Telefax: +49-30-77992-219

[www.lysoform.de](http://www.lysoform.de), e-mail: [kontakt@lysoform.de](mailto:kontakt@lysoform.de)**Lysoform Schweizerische Gesellschaft für Antisepsie AG**

Postfach 444, CH-5201 Brugg/Windisch

Telefon +41-56-441 69 81, Telefax +41-56-442 41 14

[www.lysoform.ch](http://www.lysoform.ch)

2020/08/12\_15